

Light Aluminum Premier™ Carts



Armstrong Medical

A - S M A R T® C A R T S

More Standard Features, More Flexibility,

MORE CART FOR YOUR MONEY

The Benefits Of Buying An Armstrong
ALUMINUM A-SMART® Premier™ Cart ... **Think LIGHT.**

The Benefits Of Buying **ALUMINUM** Instead of Steel

Steel

76cm Steel Cart weighs: **74.8kg**

69cm Steel Cart weighs: **72.6kg**

61cm Steel Cart weighs: **65.8kg**

Liquids can cause rust

Less maneuverability

Risk of accidental injury increases with additional weight

Aluminum

76cm Aluminum Cart weighs: **41.7kg**

69cm Aluminum Cart weighs: **38.5kg**

61cm Aluminum Cart weighs: **35.4kg**

Reduces risk of rusting

Easier maneuverability

Easier to push, turn, and **STOP**; avoid injuries to patients, visitors, nurses, and doctors

FEATURES OF A-SMART® PREMIER™ **ALUMINUM** CART

- **ALUMINUM** Carts are 45% lighter than Steel
- **Full Extension Ball Bearing Drawer Slides**
Ensures that your drawers will open and close smoothly every time.
- **Quality Aluminum Construction**
One-piece outer shell construction providing you with added strength and stability. Locking components consist of zinc plated and galvanized parts.
- **Double Side-Wall Construction**
Armstrong's A-SMART® Carts offer you the added protection of our double side-wall construction. The outer wall adds an extra layer of protection to ensure that your drawers, attached to the inner wall, will always slide smoothly.
- **Two Interchangeable Locking Systems**
- **Customize Your Locking Configuration**
Lock individual drawers, double lock drawers, use multiple types of locks on the same cart.
- **Slide-Thru Side Shelf**
- **Contour-Fit Plastic Top**
- **Four Vertical Aluminum Corner Mounting Tracks**
For easy accessory attachment (eliminates the need for holes).
- **Stabilizing Frame With Bumper**
Casters are moved out from underneath to the outside of the cart. This extra wide base provides your cart with added stability.
- **18 Colors Available**
Allows you flexibility in your carts appearance. Choose colors to specify department or function, or even use multiple colors on the same cart to differentiate drawers.
- **Highest Quality Swivel Casters (two locking, one tracking)**
Quiet, durable, German-made casters. Smoother, easier maneuvering.
- **Durable Powder Paint**
Fully powder coated paint throughout. Prevents chips, scratches, fading, and wearing.



Look for the Armstrong A-SMART® label on every cart. ***It stands for quality.***

Manufactured to ISO 9001:2008 certified standards.

Armstrong A-SMART® Push-Button Carts

Security At The Push Of A Button



Get Exactly What You Need

CREATE YOUR OWN CART



With hundreds of accessories, seventeen available colors, flexible locking options (*Push-Button, Key Locking, Breakaway Locking*), and multiple cart and drawer sizes, you can customize an A-SMART® Cart to fit any purpose or department. With so many options, you can build your own cart so that you **Get Exactly What You Need.** **A-SMART.**

Ostomy Care Cart
 Toxic Shock Cart
 Dressing Cart
 Conscious Sedation Cart
 Bronchoscopy Cart
 Eye Cart
 Block Cart
 Supply Cart
 Respiratory Cart

Trauma Cart
 Computer Cart
 Delivery Cart
 Latex-Free Cart
 Ear, Nose & Throat Cart
 Sexual Assault Exam Cart
 ER Supply Cart
 Crash Cart

Isolation Cart
 Endoscopy Cart
 Bio-Terror Cart
 Difficult Airway Cart
 Broselow Cart
 Burn Cart
 Bedside Cart
 ER Cart
 Anesthesia Cart

Phlebotomy Cart
 PICC Line Cart
 Epidural Cart
 Procedure Cart
And Many More

*Take The Virtual Tour Of Our A-SMART® Carts
 To See The Benefits of Buying An Armstrong Cart*

www.armstrongmedical.com

Table of Contents



Introduction	
<i>Standard Steel Cart Features</i>	4
<i>Premier™ All Aluminum Cart Features</i>	5
<i>Standard Steel Cart Dimensions</i>	6
<i>Premier™ All Aluminum Cart Dimensions</i>	7
PREMIUM Emergency Cart Packages	
<i>Standard Steel Emergency Cart Package</i>	8
<i>Premier™ All Aluminum Emergency Cart Package</i>	9
Emergency Accessory Packages	
<i>Standard Emergency Accessory Packages</i>	10
<i>Premier™ Emergency Accessory Packages</i>	11
Solid Color Emergency Carts	
<i>Standard Steel Solid Color Emergency Carts</i>	12
<i>Premier™ All Aluminum Solid Color Emergency Carts</i>	13
Two-Tone Emergency Carts	
<i>Standard Steel Two-Tone Emergency Carts</i>	14
<i>Premier™ All Aluminum Two-Tone Emergency Carts</i>	15
Broselow® Premier™ Aluminum Pediatric ColorCode Cart	16
Push-Button Bedside Carts	17
SUPREME Anesthesia Cart Packages	
<i>Standard Steel Anesthesia Cart Package</i>	18
<i>Premier™ All Aluminum Anesthesia Cart Package</i>	19
Anesthesia Accessory Packages	
<i>Standard Anesthesia Accessory Packages</i>	20
<i>Premier™ Anesthesia Accessory Packages</i>	21
Solid Color Mobile Workstations	
<i>Standard Steel Solid Color Key Locking Carts</i>	22
<i>Premier™ All Aluminum Solid Color Key Locking Carts</i>	23
Two-Tone Mobile Workstations	
<i>Standard Steel Two-Tone Key Locking Carts</i>	24
<i>Premier™ All Aluminum Two-Tone Key Locking Carts</i>	25
Specialty Carts and Accessories	
<i>Premier™ Locking Bars</i>	26
<i>Trellis Systems</i>	26
<i>Push-Button Locking Carts</i>	27
<i>Mini and Narrow Carts</i>	28-29
<i>ScopeSafe™ Holders & Wide Cart</i>	30
<i>Isolation Carts</i>	31
Cart Accessories	32-Inside Back Cover

Armstrong Medical Industries, Inc.

Corporate Headquarters
575 Knightsbridge Parkway • P.O. Box 700
Lincolnshire, IL 60069-0700
847/913-0101 • FAX: 847/913-0138
www.armstrongmedical.com



**Look for the Armstrong A-SMART® label
on every cart. It stands for quality.**

Manufactured to ISO 9001:2008 certified standards.

A-SMART® Setting The Standard

Breakaway Locking Cart
Breakaway Security **ABL Series**

Soft-grip Handles

Double Side-wall Construction

Ball Bearing Drawer Slides

Stabilizing Frame with Bumper prevents injuries. Casters are moved from under the cart to outside its perimeter. This extra-wide base provides your cart with added stability.

(1) Swivel
(2) Swivel w/Brake
(1) Tracking Caster

Push-Button Locking Cart
Key override comes standard on Push-Button Security **APB Series**

Key Locking Cart







APT-2 Optional Plastic Top, see page 35.



Look for the Armstrong A-SMART® label on every cart.
It stands for quality.
MADE IN THE USA

A-SMART® Standard Features:

- Ball bearing drawer slides
- High-quality swivel casters (two locking, one tracking)
-  Casters
- Double Side-wall construction (see right)
- Soft-grip handles
- Stabilizing frame with bumper
- 18 Colors
- Durable powder paint
- All carts are latex free

Double Side-Wall Construction

Armstrong's A-SMART® Carts offer you the added protection of our double side-wall construction. The outer wall adds an extra layer of protection to ensure that your drawers, attached to the inner wall, will always slide smoothly.

Take The Strain Off Your Back. Think **LIGHT.**



Customize your locking configuration; lock individual drawers, double-lock drawers, use multiple types of locks on the same cart.

Key Locking Bar

Has all the same features as the breakaway locking cart. Security Locking Bar PKL Series



Look for the Armstrong A-SMART® label on every cart.

It stands for quality.

MADE IN THE USA

See the Advantages of A-SMART® PREMIER™ Aluminum Carts

- Lightweight, aluminum construction
- Two interchangeable locking systems
- Four vertical Aluminum corner mounting tracks for easy accessory attachment
- Slide-thru side shelf
- Contour-fit plastic top
- High-quality swivel casters (two locking, one tracking)
- Stabilizing frame with bumper
- Double side-wall construction
- Soft-grip handles
- 18 colors
- Durable powder paint
- Ball bearing drawer slides
- All carts are latex free



Double Side-Wall Construction

Armstrong's A-SMART® Carts offer you the added protection of our double side-wall construction. The outer wall adds an extra layer of protection to ensure that your drawers, attached to the inner wall, will always slide smoothly.

Standard Steel Cart Features

How to Customize An A-SMART® Standard Steel Cart:

1. CHOOSE FRAME SIZE

What kind of locking system do you need? How much drawer space do you want?

ABL Series: Breakaway Locks



ABL-6

6 Drws. } = 76cm drawer space

ABL-6 = 6 Drws. with 76cm dr. space (shown above)
ABL-4 = 4 Drws. with 61cm dr. space
ABL-5 = 5 Drws. with 61cm dr. space


6-Dr. ABL Cart Dim.:

Overall: 117.5cm High x 63.5cm Deep x 85.1cm Wide (includes bumper/handles/casters).
Drawer size: 43.3cm Deep x 56.5cm Wide
6-Drawer: 7.6cm (4), 15.6cm (1), 31.8cm (1).
Weight: 74.8kg

4-Dr. and 5-Dr. ABL Cart Dim.

Overall: 101cm High x 65.3cm Deep x 85.1cm Wide (includes bumper/handles/casters).
Drawer size: 43.3cm Deep x 56.5cm Wide
4-Drawer: 7.3cm (2), 15.6cm (1), 32.2cm (1).
5-Drawer: 7.3cm (3), 15.6cm (1), 23.8cm (1).
Weight: 65.8kg

AKL Series: Key Locks



AKL-6

6 Drws. } = 76cm drawer space

AKL-6 = 6 Drws. with 76cm dr. space (shown above)
AKL-4 = 4 Drws. with 76cm dr. space
AKL-5 = 5 Drws. with 61cm dr. space

4-Dr. and 6-Dr. AKL Cart Dim.:


Overall: 117.5cm High x 63.5cm Deep x 85.1cm Wide (includes bumper/handles/casters).
Drawer size: 43.3cm Deep x 56.5cm Wide
4-Drawer: 7.3cm (1), 23.8cm (3).
6-Drawer: 7.3cm (3), 15.6cm (2), 23.8cm (1).
Weight: 74.8kg

5-Dr. AKL Cart Dim.:

Overall: 101cm High x 65.3cm Deep x 85.1cm Wide (includes bumper/handles/casters).
Drawer size: 43.3cm Deep x 56.5cm Wide
5-Drawer: 7.3cm (3), 15.6cm (1), 23.8cm (1).
Weight: 65.8kg

APB Series:

Push-Button Locks w/Key Override



APB-6

6 Drws. } = 76cm drawer space

APB-6 = 6 Drws. with 76cm dr. space (shown above)
APB-5 = 5 Drws. with 61cm dr. space

Standard Push-Button Locking Carts, page 27


6-Dr. APB Cart Dim.:

Overall: 117.5cm High x 63.5cm Deep x 85.1cm Wide (includes bumper/handles/casters).
Drawer size: 43.3cm Deep x 56.5cm Wide
6-Drawer: 7.3cm (3), 15.6cm (2), 23.8cm (1).
Weight: 74.8kg

5-Dr. APB Cart Dim.:

Overall: 101cm High x 65.3cm Deep x 85.1cm Wide (includes bumper/handles/casters).
Drawer size: 43.3cm Deep x 56.5cm Wide
5-Drawer: 7.3cm (3), 15.6cm (1), 23.8cm (1).
Weight: 65.8kg

18 Designer Colors



2. Pick Your Options


Choose your drawer fronts...

Choose from any one of our fourteen designer colors on the right. If you want the ABL-4 as shown with Teal drawer fronts, your cart would be ABL-T-4.

Four drawer sizes available:

7.3cm 15.6cm 23.8cm 31.8cm
 If standard drawer arrangements don't suit your needs, slides can be moved 7.3cm up or down (additional charge per cart).

18 Designer Colors



How To Customize An A-SMART® All Aluminum Premier™ Cart:

1. Choose Frame Size




61cm (5-Drw) CART DIMENSIONS:

Overall: 99.1cm High x 63.5cm Deep x 83.2cm Wide (includes bumper, handles, and casters).

Drawer size: 43.3cm Deep x 56.5cm Wide.

5-Drawer: 7.3cm (3), 15.6cm (1), 23.8cm (1).

Weight: 35.4kg



68.6cm (6-Drw) CART DIMENSIONS:

Overall: 107.3cm High x 63.5cm Deep x 83.2cm Wide (includes bumper, handles, and casters).

Drawer size: 43.3cm Deep x 56.5cm Wide.

6-Drawer: 7.3cm (4), 15.6cm (1), 23.8cm (1).

Weight: 38.6kg



76.2cm (6-Drw) CART DIMENSIONS:

Overall: 114.9cm High x 63.5cm Deep x 83.2cm Wide (includes bumper, handles, and casters).

Drawer size: 43.3cm Deep x 56.5cm Wide.

6-Drawer: 7.3cm (3), 15.6cm (2), 23.8cm (1).


Weight: 41.7kg

2. Choose a Lock or Have Both

PKL Series: Key Locks



PBL Series: Breakaway Locks



Locks fit either or both corner glides



3. Pick Your Options

Choose your drawer fronts...

Choose from any one of the eighteen designer colors on the left. To get the PKL-24 as shown with **Hunter Green** drawer fronts, order the PKL-HG-24. To get the PBL-27 as shown with **Red** drawer fronts, order the PBL-AR-27.

Four drawer sizes available:

7.3cm 15.6cm 23.8cm 31.8cm

If standard drawer arrangements don't suit your needs, slides can be moved 7.3cm up or down (additional charge per cart).



Look for the Armstrong A-SMART® label on every cart. *It stands for quality.*

MADE IN THE USA

**ACCESSORIES,
Pages 32-Inside Back Cover**

Standard Premium Emergency Cart Package

BREAKAWAY LOCKING

Standard Premium Emergency Package AR-6 with APE-1X Emergency Package

**Double
Side-Wall
Construction**

AUH-1
Utility Hooks, 2/pkg.
(not pictured)

APT-2
Plastic Top

ATW-3
Soft Grip
Handles

ADT-1
Drug Tray with Lid

ADV-10
Modular Dividers
for ADT-1

AC-670S (4)
7.3cm High Drug
Trays with Dividers

AC-672S (2)
12.7cm High Drug
Trays with Dividers

**Quiet, Easy-roll,
Ball Bearing
Drawer Slides**



AOS-1
Hospital-grade
Electrical Outlet Strip

AOB-1
Oxygen Tank Bracket

PIV-1
I.V. Pole

ADS-1N
Defib Swivel
Shelving Unit

ASN-3
Collapsible Side
Shelf

ALH-12
Label Holders

APS-5
Plastic Seals,
100/pkg.
(not pictured)

**All Drawers
Have Full
Extension**

ASU-1
Suction Unit
Holder

ACB-1
Cardiac Board
with Brackets
(not pictured)


ASF-2
Standard
Stabilizing Frame
with Bumpers

See page 16 for the
Broselow® ColorCode™ Cart

NOTE: Carts require some assembly.

Premier™ Premium Emergency Package PBL-AR-30 with PAPE-1X Emergency Package

Double Side-Wall Construction



PCB-1
Cardiac Board w/Brackets

See page 26
for Premier™
Locking Options.

Standard Emergency Cart Accessory Packages

BREAKAWAY LOCKING

Deluxe Emergency Accessory Package

EB-1 Deluxe Accessory Package *(Does not include cart)*

Complete Package as shown, EB-1 with AR-6 Six-Drawer 76cm Emergency Cart

- ASN-3 Collapsible Side Shelf
- PIV-1 I.V. Pole with Brackets
- APT-2 Plastic Top
- PSR-1 Heavy-duty Railing
- AOB-1 Oxygen Tank Bracket
- ACB-1 Cardiac Board with Brackets
- APS-5 Plastic Breakaway Seals, 100/pkg.
- AUH-1 Utility Hooks, 2/pkg
- Two AC-672S 12.7cm High Trays with four AC-685S (short) Dividers & two AC-690S (long) Dividers (fills one drawer)
- Four AC-670S 7.3cm High Trays with eight AC-675S (short) Dividers and four AC-680S (long) Dividers (fills two drawers)



Modular Drawer Dividers
see page 40



Look for the
Armstrong A-SMART® label
on every cart.
It stands for quality.


MADE IN THE USA

Economy Emergency Accessory Package

BE-1 Economy Accessory Package *(Does not include cart)*

Complete Package as shown, BE-1 with ABL-AR-6 Six-Drawer 76cm Emergency Cart

- PIV-1 I.V. Pole with Bracket
- PSR-1 Heavy-duty Railing
- ACB-1 Cardiac Board with Brackets
- APT-2 Plastic Top
- APS-5 Plastic Breakaway Seals, 100/pkg.
- AOB-1 Oxygen Tank Bracket



Tilt Bins, see page 37



NOTE: Carts require some assembly.

See page 4 for Standard Cart features
and page 6 for dimensions.


BREAKAWAY LOCKING

Deluxe Emergency Accessory Package

PEB-1 Deluxe Accessory Package *(Does not include cart)*

Complete Package as shown, PEB-1 with PAR-30 Six-Drawer 76cm All Aluminum Emergency Cart

- PIV-1 I.V. Pole with Brackets
- PSR-1 Tubular Steel Railing
- AOB-1 Oxygen Tank Bracket
- PCB-1 Cardiac Board with Brackets
- APS-5 Breakaway Seals, 100/pkg.
- Two AC-672S 12.7cm High Trays with four AC-685S (short) Dividers and two AC-690S (long) Dividers
- Four AC-670S 7.3cm High Trays with eight AC-675S (short) Dividers and four AC-680S (long) Dividers



ADV-1 Modular Drawer Divider Kit, see page 39

Mount either of these accessory packages on any A-SMART® Premier™ All Aluminum Emergency Cart



Look for the Armstrong A-SMART® label on every cart. *It stands for quality.*

MADE IN THE USA

Economy Emergency Accessory Package

PBE-1 Economy Accessory Package *(Does not include cart)*

Complete Package as shown, PBE-1 with PBL-AB-30 Six-Drawer 76cm All Aluminum Emergency Cart

- PIV-1 I.V. Pole with Brackets
- PSR-1 Tubular Steel Railing
- AOB-1 Oxygen Tank Bracket
- PCB-1 Cardiac Board with Brackets *(shown above)*
- APS-5 Breakaway Seals, 100/pkg.



See page 5 for Premier™ Cart features and page 7 for dimensions

Standard Steel Solid Color Emergency Carts

BREAKAWAY LOCKING

Manufactured To ISO 9001:2008 Certified Standards.

A-SMART®.



ABL-B-6

Six-Drawer 76cm Standard Emergency Cart

ABL-B-6 Emergency Cart, Beige
Drawer sizes: 7.3cm (4), 15.6cm, 31.8cm

Double-wall construction and ball bearing drawer slides add to the quality of all Armstrong Carts.



ABL-B-5

Five-Drawer 61cm Standard Emergency Cart

ABL-B-5 Emergency Cart, Beige
Drawer sizes: 7.3cm (3), 15.6cm, 23.8cm

Built-in, stabilizing frame/bumper gives Armstrong Carts the stability for emergency situations.



ABL-B-4

Four-Drawer 61cm Standard Emergency Cart

ABL-B-4 Emergency Cart, Beige
AR-4 Emergency Cart, Red

(not pictured)

Drawer sizes: 7.3cm (2), 15.6cm, 31.8cm



AR-6

Six-Drawer 76cm Standard Emergency Cart

AR-6 Emergency Cart, Red
Drawer sizes: 7.3cm (4), 15.6cm, 31.8cm

Double-wall construction and ball bearing drawer slides add to the quality of all Armstrong Carts.



AR-4

Four-Drawer 61cm Standard Emergency Cart

AR-4 Emergency Cart, Red
Drawer sizes: 7.3cm (2), 15.6cm, 31.8cm



AR-27

Six-Drawer 69cm Standard Emergency Cart

AR-27 Emergency Cart, Red
Overall Dim.: 109.2cm High x 63.5cm Deep x 85.1cm Wide
(includes bumper, handles, and casters)
Drawer sizes: 7.3cm (4), 15.6cm, 23.8cm
Weight: 72.6kg

See page 4 for Standard Cart features and page 6 for dimensions.

NOTE: Carts require some assembly.

Drawers Always Slide Smoothly With Our Double Side-Wall Construction. . . **Think LIGHT.**



PBL-B-27

Six-Drawer 69cm Premier Emergency Cart

PBL-B-27 Emergency Cart, Beige
Drawer sizes: 7.3cm (4), 15.6cm, 23.8cm



PBL-B-24

Five-Drawer 61cm Premier Emergency Cart

PBL-B-24 Emergency Cart, Beige
Drawer sizes: 7.3cm (3), 15.6cm, 23.8cm



PAR-30

Six-Drawer 76cm Premier Emergency Carts

PAR-30 Emergency Cart, Red
PBL-B-30 Emergency Cart, Beige
(not pictured)
Drawer sizes: 7.3cm (3), 15.6cm (2), 23.8cm



AC-670S and AC-672S,
see page 40




PSR-1 Tubular Railing,
see page 35



Look for the
Armstrong A-SMART® label
on every cart.
It stands for quality.

MADE IN THE USA

Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster  Casters
- Tracking caster
- Breakaway locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks

Standard Steel Two-Tone Emergency Carts

BREAKAWAY LOCKING

Armstrong A-SMART® Two-Tone Carts
come with your choice of 18 colored drawer fronts.



ABL-L-5

**Five-Drawer 61cm
Standard Emergency Cart**

ABL-L-5 Emergency Cart,
Lavender
Drawer sizes: 7.3cm (3), 15.6cm,
23.8cm



ABL-T-4

**Four-Drawer 61cm
Standard Emergency Cart**

ABL-T-4 Emergency Cart,
Teal
Drawer sizes: 7.3cm (2), 15.6cm,
31.8cm



ABL-TC-6

**Six-Drawer 76cm Standard
Emergency Cart**

ABL-TC-6 Emergency Cart,
Terra Cotta
Drawer sizes: 7.3cm (4), 15.6cm,
31.8cm



Look for the Armstrong
A-SMART® label on every
cart. *It stands for quality.*

MADE IN THE USA

**Three-Drawer 61cm
Standard Emergency Cart**

ABL-B-3 Emergency Cart
Beige

Drawers: 15.6cm, 23.8cm (2)



ABL-B-3

18 Designer Colors

SG - Sage Green

MG - Moss Gray

TC - Terra Cotta

SB - Slate Blue

L - Lavender

M - Mauve

Y - Yellow

CB - Taupe

T - Teal

B - Beige

AB - Crash Cart Blue

CR - Cranberry

BG - Blue Grey

HG - Hunter Green

AA - Deep Blue

AR - Red

S - Silver

CW - Cherry Woodgrain



**AHP-3 Panel for Cart-mounted
DUET Suction Unit**
(Retention Bracket not included),
see page 34



PSR-1 Tubular Railing,
see page 35




ASN-3 Standard Side Shelf,
see page 35

Carts Built To Last Longer!

Think LIGHT.

18 Designer Colors



PBL-AR-30

Six-Drawer 76cm Premier Emergency Cart

PBL-AR-30 Emergency Cart, Red
Drawer sizes: 7.3cm (3), 15.6cm (2), 23.8cm



PBL-SG-27

Six-Drawer 69cm Premier Emergency Cart

PBL-SG-27 Emergency Cart, Sage Green
Drawer sizes: 7.3cm (4), 15.6cm, 23.8cm



PBL-SB-24

Five-Drawer 61cm Premier Emergency Cart

PBL-SB-24 Emergency Cart, Slate Blue
Drawer sizes: 7.3cm (3), 15.6cm, 23.8cm



AMT-2N Modular Tray, see page 38




PDS-1 One shelf Unit, Deep, see page 36



Look for the Armstrong A-SMART® label on every cart. It stands for quality.

MADE IN THE USA

Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster  Casters
- Tracking caster
- Breakaway locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks

Premier™ Aluminum ColorCode™ Cart

BREAKAWAY LOCKING

Not MRI Compatible

Use The ColorCode™ Cart In Pediatric Resuscitation Emergencies.



PBL-PC-9


Use OPTIONAL 7.3cm drawer locks so that only the accessed drawers need restocking.



TENTE Casters

PBL-PC-9 BASIC A-SMART®
Premier™ Broselow® ColorCode™
Cart (pictured above)

Includes Nine-Drawer A-SMART Premier Breakaway Lock Cart (76cm total drawer space), slide-thru side shelf, 76cm locking bar, stabilizing frame/bumper system, seamless plastic top, soft-grip push handles, three accessory attachment panels, and all-swivel casters (two locking and one tracking).



PBL-PC-9 shown above with OPTIONAL Accessory Package
PC-9E: Includes AOB-1 Oxygen Tank Bracket, ASU-1 Suction Unit Holder, PIV-1


IV Pole, PCB-1 Cardiac Board (not shown), PBR-3 (8) 7.3cm Breakaway Locking Bar, PBR-6 (1) 15.6cm Breakaway Locking Bar, and APS-5 Plastic Seals.

Greater Continuity, Accuracy, and Efficiency

PBR-3 7.3cm Breakaway Locking Bar
PBR-6 15.6cm Breakaway Locking Bar

Alternative locking options help simplify inventory or limit access to sensitive materials. See page 26 for more locking options.

18 Designer Colors



NURSE SERVER/BEDSIDE Carts



Aluminum Push-Button Locking Narrow Cart

AMC-3Color Six-Drawer 76cm Push-Button Locking Narrow Cart

Beige shell with eleven drawer colors available (see page 29).

Features:

- Aluminum construction; weighs only 38.5kg
- Stabilizing frame w/ bumper
- Plastic top
- Four 12.7cm swivel casters (two with brake and one tracking)
- Push-button locking with key override

Dimensions:

Overall: 111.8cm High x 62.9cm Deep x 62.2cm Wide (includes bumper and casters).

Drawer configuration: 76cm of total drawer height.

Drawer size: 43.2cm Deep x 38.1cm Wide.

Drawer heights: 7.3cm (3), 15.6cm (2), 23.8cm (1).

Aluminum Five-Drawer 61cm Push-Button Locking Cart

PPB-Color-5 Five-Drawer 61cm Push-Button Locking Cart

Beige shell with eighteen drawer colors available.

Drawer sizes: 7.3cm (3), 15.6cm (1), 23.8cm (1).

Weight: 68kg

Four 12.7cm swivel casters (two with brake and one tracking)

See page 27 for more details.

Steel Five-Drawer 61cm Push-Button Locking Mini Cart

AMC-1Color-PB Five-Drawer 61cm Push-Button Locking Mini Cart with Key Override

Beige shell with eleven drawer colors available (see page 28).

Overall: 88.9cm High x 52.1cm Deep x 53.3cm Wide.

Drawer size: 43.2cm Deep x 38.1cm Wide.

Drawer heights: 7.3cm (3), 15.6cm (1), 23.8cm (1).

Weight: 46.7kg

Four 7.6cm swivel casters (two with brake and one tracking)

Steel Four-Drawer 53cm Push-Button Locking Mini Cart

AMC-21Color-PB Four-Drawer 53cm Push-Button Locking Mini Cart with Key Override

Beige shell with eleven drawer colors available (see page 28).

Overall: 80cm High x 52.1cm Deep x 53.3cm Wide.

Drawer size: 43.2cm Deep x 38.1cm Wide.

Drawer heights: 7.3cm (2), 15.6cm (1), 23.8cm (1).


Four 7.3cm swivel casters (two with brake and one tracking)

Standard Steel Supreme Anesthesia Cart Package

KEY-LOCKING

STANDARD SUPREME ANESTHESIA PACKAGE

Complete Steel Cart with entire Accessory Package as shown
AKL-AA-6 with ASA-1X Package



Look for the
Armstrong A-SMART® label
on every cart.
It stands for quality.


NOTE: Carts require some assembly.

Premier™ Aluminum Supreme Anesthesia Cart Package

KEY-LOCKING
Not MRI Compatible

PREMIER™ SUPREME ANESTHESIA PACKAGE

Complete Aluminum Cart with entire Accessory Package as shown
PKL-AA-30 with PASA-1X Package



See page 26 for Premier™ Locking Options.

Look for the
Armstrong A-SMART® label
on every cart.
It stands for quality.

MADE IN THE USA

NOTE: Carts require some assembly.

These Accessory Packages Mount To Any Standard Steel Cart

Deluxe Anesthesia Accessory Package

DAB-1 Deluxe Accessory Package (Does not include cart)

Complete Package as shown, DAB-1 and APB-AA-6 Six-Drawer 76cm Push-Button Locking Cart

- ATR-6N Trellis System
- ATB-4X Four-bin Tilt Bin Organizer
- ATB-5X Five-bin Tilt Bin Organizer
- PSR-1 Tubular Railing
- ASN-3 Collapsible Side Shelf
- AWC-1 Waste Container
- AAB-1 Accessory Bar
- APT-2 Plastic Top



PIV-1 I.V. Pole with Brackets, see page 35



Look for the Armstrong A-SMART® label on every cart. *It stands for quality.*

MADE IN THE USA



APB-AA-6

Push-Button Carts, see page 27



APB-B-6

Economy Anesthesia Accessory Package

SP-1B Economy Accessory Package (Does not include cart)

Complete Package as shown, SP-1B and APB-B-6 Six-Drawer 76cm Push-Button Locking Cart

- APT-2 Plastic Top
- AUH-1 Utility Hooks, 2/pkg.
- PSR-1 Tubular Railing
- Two AC-672S 12.7cm Drug Trays with four Short Dividers and two Long Dividers (fills one drawer)
- Four AC-670S 7.3cm Drug Trays with eight Short Dividers and four Long Dividers (fills two drawers)



Modular Drawer Trays, see page 38

NOTE: Carts require some assembly.

Accessory Packages To Lighten The Load. . .

Think LIGHT.

Deluxe Anesthesia Accessory Package

PDAB-1 Premier Deluxe Accessory Package *(Does not include cart)*

Complete Package as shown, PDAB-1 and PKL-SB-30 Six-Drawer 76cm All Aluminum Mobile Workstation

- PNS-1 Shelf Unit
- ATB-5X Five-bin Tilt Bin Organizer
- AGB-1 Glove Box
- PSB-1 Sharps Container Bracket
- AWC-1 Waste Container



PKL-SB-30



ScopeSafe™
*Protect Your
Delicate
Scopes*

ScopeSafe™
Locking Scope
Holder,
see page 30

Economy Anesthesia Accessory Package

PSP-1B Premier Economy Accessory Package
(Does not include cart)

Complete Package as shown, PSP-1B and PKL-CW-30
Six-Drawer 76cm All Aluminum Mobile Workstation

- AGB-1 Glove Box
- AWC-1 Waste Container
- Four AC-670S 7.3 cm Drug Trays with eight Short Dividers and four Long Dividers
- Two AC-672S 12.7 cm Drug Trays with four Short Dividers and two Long Dividers



Premier™
Locking
Bars,
see page 26



PKL-CW-30

See page 5 for Premier™ Cart features and page 7 for dimensions.

Standard Steel Solid Color Mobile Workstations

KEY-LOCKING

Solid Colors, Solid Features

All A-SMART® Key Locking Carts come with double side-wall construction, ball bearing drawer slides, a standard stabilizing frame/bumper, and four all-swivel casters (two locking and one tracking).



Six-Drawer 76cm Standard Mobile Workstations

AA-6 Mobile Workstation,
Deep Blue

Drawer sizes: 7.3cm (3), 15.6cm (2),
23.8cm



Six-Drawer 76cm Standard Mobile Workstations

AN-6 Mobile Workstation,
Beige

Drawer sizes: 7.3cm (3), 15.6cm (2),
23.8cm



Five-Drawer 61cm Standard Mobile Workstations

AN-5 Mobile Workstation,
Beige

Drawer sizes: 7.3cm (3), 15.6cm,
23.8cm



Look for the
Armstrong A-SMART® label
on every cart.

It stands for quality.

Want more accessories? See pages 32-Inside Back Cover



PIV-1 I.V. Pole with
Brackets, see page 35



ATR-6N Standard Trellis
Rail System, see page 26



ASN-3 Side Shelf,
see page 35



APT-2 Standard Plastic
Top, see page 35

Standard Anesthesia Packages, pages 18, 20
Cart Accessories, pages 32-Inside Back Cover

See page 4 for Standard Cart features and
page 6 for dimensions.

Premier™ Carts Are Easier To Push and STOP. . .

Think LIGHT.



PKL-B-27

Six-Drawer 69cm Premier Mobile Workstation

PKL-B-27 Mobile Workstation,
Beige

Drawer sizes: 7.3cm (4), 15.6cm, 23.8cm



AEF-6 Plastic Divider Set
for 6" drawer, see page 38



PKL-B-24

Five-Drawer 61cm Premier Mobile Workstation

PKL-B-24 Mobile Workstation,
Beige

Drawer sizes: 7.3cm (3), 15.6cm,
23.8cm



PDS-1 One-shelf Unit, Deep,
see page 36



PKL-B-30

Six-Drawer 76cm Premier Mobile Workstation

PKL-B-30 Mobile Workstation,
Beige


Drawer sizes: 7.3cm (3), 15.6cm (2),
23.8cm



Look for the
Armstrong A-SMART® label
on every cart.
It stands for quality.

Visit Us At:
[www.armstrongmedical.11`](http://www.armstrongmedical.11)

Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster  Casters
- Tracking caster
- Breakaway locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks


Standard Steel Two-Tone Mobile Workstations

KEY-LOCKING

Colors To Match Your Decor

When your colors are added to the quality construction of an Armstrong A-SMART® Cart, you've got a Mobile Workstation like no other. **A-SMART®**

18 Designer Colors



AKL-BG-4

Four-Drawer 76cm Standard Mobile Workstation

AKL-BG-4 Mobile Workstation, Blue Grey

Drawer sizes: 7.3cm, 23.8cm (3)



AKL-SB-6

Six-Drawer 76cm Standard Mobile Workstation

AKL-SB-6 Mobile Workstation, Slate Blue

Drawer sizes: 7.3cm (3), 15.6cm (2), 23.8cm



AKL-MG-5

Five-Drawer 61cm Standard Mobile Workstation

AKL-MG-5 Mobile Workstation, Moss Gray

Drawer sizes: 7.3cm (3), 15.6cm, 23.8cm



Look for the Armstrong A-SMART® label on every cart. **It stands for quality.**

MADE IN THE USA



Standard Side Shelf, see page 35



Modular Drawer Trays, see page 38



Wire Baskets, see page 34


Standard Anesthesia Packages, pages 18, 20

Accessories, pages 32-Inside Back Cover

Maneuverability At It's Best. . .

Think LIGHT.

18 Designer Colors



Six-Drawer 69cm Premier Mobile Workstation

PKL-CB-27 Mobile Workstation, Taupe

Drawer sizes: 7.3cm (4), 15.6cm, 23.8cm

Six-Drawer 76cm Premier Mobile Workstation

PKL-HG-30 Mobile Workstation, Hunter Green

Drawer sizes: 7.3cm (3), 15.6cm (2), 23.8cm

Four-Drawer 76cm Premier Mobile Workstation

PKL-Y-4 Mobile Workstation, Yellow

Drawer sizes: 7.3cm, 23.8cm (3)



Look for the Armstrong A-SMART® label on every cart.

It stands for quality.

MADE IN THE USA

www.armstrongmedical.com


Five-Drawer 61cm Premier Mobile Workstation

PKL-SB-24 Mobile Workstation, Slate Blue

Drawer sizes: 7.3cm (3), 15.6cm, 23.8cm



Aluminum Cart Features:

- Ball bearing drawer slides
- Slide-thru side shelf
- Plastic top
- Soft-grip handles (2)
- Accessory panels (3)
- Stabilizing frame with bumper
- Locking casters (2)
- Swivel caster  Casters
- Tracking caster
- Key locking bar (length of bar is equal to total drawer space)
- Aluminum mounting tracks

Premier™ Locking Bars and Trellis Systems

Restock Only The Accessed Drawers. Save Money. Save Time.

Don't waste money checking crash cart drawers that haven't even been opened. With A-SMART® Premier™ Carts, individual drawer seals can be arranged so that only the accessed drawers need to be checked for inventory, a small investment that pays off with every use.

- Double lock medication drawers
- Improve equipment security
- Ready for changing needs

When carts are shared by departments or personnel, A-SMART Premier's locking options allow each to be locked with a separate key. Never wonder about who else has access to your equipment.

Maybe you don't need this kind of flexibility now, but you might later. With A-SMART Premier, you can update your locking system at any time. Need to change a mobile workstation into a crash cart? No problem. Need to add a key lock to an existing cart? Done. Only A-SMART Premier offers you that kind of flexibility.



Locks may be placed on either corner glide for flexible security. Double locking drawers can help control inventory and access.



Key Locking Bars

Different sized locking bars help control access.


Size	Prod#
7.3cm	PKE-3
15.6cm	PKE-6
23.8cm	PKE-9
61cm	PKE-24
69cm	PKE-27
76cm	PKE-30



Breakaway Locking Bars

Mix lock sizes and types for precise inventory control.

Size	Prod#
7.3cm	PBR-3
15.6cm	PBR-6
23.8cm	PBR-9
38.1cm	PBR-15
53.3cm	PBR-21
61cm	PBR-24
69cm	PBR-27
76cm	PBR-30




PTR-1 PREMIER™ Trellis System
Adds up to 53.3cm to cart height. Bars are height-adjustable. Includes two uprights and three cross bars.

PTR-5 PREMIER™ Trellis System with Hooks (not pictured)
Adds up to 53.3cm to cart height. Bars are height-adjustable. Includes two uprights and three cross bars. Includes one bar with five hooks in the top rail.



ATR-5N STANDARD Trellis Rail System with Five Hooks
Same features as the ATR-6N (see right), but with five hooks on the top rail for hanging supplies.

PHB-1 Additional Hanging Bar for Standard or Premier Trellis Systems




ATR-6N STANDARD Trellis Rail System
The top bar is 54.6cm above top of cart.

ATR-9N Wide Cart Trellis Rail System
Fits Armstrong's Wide AKW Series Carts (see page 30).

NEW Aluminum Electronic Touchpad and Steel Push-Button Carts

TOUCHPAD AND PUSH-BUTTON LOCKING
Not MRI Compatible

18 Designer Colors



NEW Aluminum ELECTRONIC Touchpad Cart

PPE-Color-76cm Six-Drawer Aluminum Electronic Touchpad Cart
Drawer sizes: 7.3cm (3), 15.6cm (2), 23.8cm (1)
Weight: 41.7 kg.

PPE-Color-61cm Five-Drawer Aluminum Electronic Touchpad Cart
Drawer sizes: 7.3cm (3), 15.6cm, 23.8cm. (1)
Weight: 34.4 kg.

Features:

- Simple touchpad entry with key override
- 1-8 digits
- Operates on alkaline batteries
- User code is stored when battery power is lost or removed
- Audible low battery indicator
- Continuous alarm sounds when code is incorrectly entered four times
- Plastic top
- Accessory panels (3)
- Stabilizing frame with bumper
- Aluminum mounting tracks
- 18 colors available

NOTE: No auto-locking features (control knob must be turned manually to lock and unlock). Carts require some assembly.

* Please note for the Armstrong Electronic Touchpad Carts, we suggest using alcohol only for the electronic keypad surface. Please completely dry keypad and surrounding area after alcohol application.



PPE-AA-30

PPE-AA-30 shown with optional PSR-1 Railing, ATD-2 Tape Dispenser, PTR-1 Trellis System, and ATB-4X & ATB-5X Tilt Bins



APB-SB-6

Standard Steel Push-Button Locking Carts

APB-Color-6 Standard Six-Drawer 76cm Push-Button Locking Cart
Drawer sizes: 7.3cm (3), 15.6cm (2), 23.8cm.
Weight: 77kg

APB-Color-5 Standard Five-Drawer 61cm Push-Button Locking Cart
Drawer sizes: 7.3cm (3), 15.6cm, 23.8cm.
Weight: 68kg

This system comes with a key override to bypass Push-Button System if needed.

18 colors available.

APB-SB-6 shown at left with optional APT-2 Plastic Top and PSR-1 Tubular Railing.

Steel Mini Carts

PUSH-BUTTON LOCKING & KEY LOCKING

Spacesaver

Steel Two-Tone Mini Carts

Ideal For Treatments, Procedures, & the Bedside. *Small enough to put under counters.*

Carts below have a beige shell with eleven drawer colors available:

Compact, five-drawer Mini Carts include a plastic top, ball bearing drawer slides, corner bumper guards, and rugged 8.9cm swivel casters (two locking).

Mini Carts also features colored drawer fronts in a beige frame.

Overall: 88.9cm High x 52.1cm Deep x 53.3cm Wide.

Drawer size: 43.2cm Deep x 38.1cm Wide.

Drawer heights: 7.3cm (3), 15.6cm, 23.8cm

Weight: 46.7kg

Corner bumpers do not come assembled.

See page 29 for Mini Cart accessories.



AMC-1SG-PB



AMC-21CW-PB



AMC-1MG

11 Designer Colors

SG - Sage Green

MG - Moss Grey

TC - Terra Cotta

AA - Deep Blue

B - Beige

BG - Blue Grey

HG - Hunter Green

T - Teal

Y - Yellow

CB - Taupe

CW - Cherry Woodgrain

Steel Push-Button Locking Mini Cart

AMC-1Color-PB Five-Drawer 61cm Push-Button Locking Mini Cart

www.armstrongmedical.com

Steel Push-Button 53cm Locking Mini Cart

AMC-21Color-PB Four-Drawer 53cm Push-Button Locking Mini Cart

Overall: 80cm High x 52.1cm Deep x 53.3cm Wide.

Drawer size: 43.2cm Deep x 38.1cm Wide.


Drawer heights: 7.3cm (2), 15.6cm, 23.8cm.

Weight: 38.5kg

Steel Key Locking Mini Cart

AMC-1Color Five-Drawer 61cm Key Locking Mini Cart

11 Designer Colors



Great For All Those Tight Spaces

Features:

- Lightweight aluminum construction; weighs only 38.6kg
- Stabilizing frame with bumper
- Four 12.7cm swivel casters (two with brake and one tracking)
- Plastic top
- Eleven colors available

Dimensions:

Overall: 111.8cm High x 62.9cm Deep x 62.2cm Wide (includes bumper and casters).

Drawer configuration: 76cm of total drawer height.

Drawer size: 43.2cm Deep x 38.1cm Wide.

Drawer heights: 7.3cm (3), 15.6cm (2), 23.8cm (1).

See below for Narrow Cart accessories.

**Glove Dispensers,
Inside Back Cover**



AMC-2HG

Aluminum 76cm Key Locking Narrow Cart

AMC-2Color Six-Drawer 76cm Key Locking Narrow Cart



AMC-3BG

Aluminum 76cm Push-Button Locking Narrow Cart

AMC-3Color Six-Drawer 76cm Push-Button Locking Narrow Cart

Optional **MINI** & **NARROW** Cart Accessories

APH-1 Push Handle

MDM-1B Keyboard Drawer Weight capacity: 13.6kg

ATR-7N Three-rail Trellis System

ATR-8N Three-rail Trellis System with Four Hooks

AMT-4N Full Drawer Tray with Nine Ampule Dividers, Ampule Tray, Rail, and

Four Adjustable Dividers (Fits one per drawer), see page 38.

AMT-5N Full Drawer Tray with Two Rails and Four Adjustable Dividers (Fits one per drawer), see page 38.

AEF-9M Plastic Divider Set for 23.8cm

Mini Cart Drawer

MDM-1B Keyboard Drawer

AEF-6M Plastic Divider Set for 15.6cm Mini Cart Drawer

AEF-3M Plastic Divider Set for 7.3cm Mini Cart Drawer

Also compatible with all accessories that fit on the sides of our Standard Mobile Workstations. See pages 32-Inside Back Cover for more accessories.

Aluminum Wide Cart and ScopeSafe™

KEY-LOCKING



AKE-B-6

NEW Aluminum ELECTRONIC Touchpad WIDE Cart

AKE-B-6 6-Drawer 76cm Electronic Touchpad Wide Cart, Beige

A-SMART® Wide Carts feature:

- Simple touchpad entry with key override
- 1-8 digits
- Operates on alkaline batteries
- User code is stored when battery power is lost or removed
- Audible low battery indicator
- Continuous alarm sounds when code is incorrectly entered four times
- Extra-wide drawers: 7" wider than standard cart drawers
- Seamless, one-piece plastic top
- Standard stabilizing frame with bumper
- 5" all-swivel casters (two locking, one tracking)
- Beige drawer fronts in beige frame
- Carts weigh 53.5 kg.

Dimensions:

Overall size: 115.6cm H x 64.5cm D x 105.8cm W (includes bumper, handles, and casters).

Actual Top Dimensions: 82.5cm W x 43.2cm D

Available drawer configuration:

Six-Drawer: 7.3cm (3), 15.6cm (2), 23.8cm (1).

NOTE: No auto-locking features (control knob must be turned manually to lock and unlock). Carts require some assembly.

A-SMART WIDE CART SYSTEMS ACCEPT ALL THE ACCESSORIES THAT FIT ON THE SIDES OF OUR STANDARD MOBILE WORKSTATIONS.

** Please note for the Armstrong Electronic Touchpad Carts, we suggest using alcohol only for the electronic keypad surface. Please completely dry keypad and surrounding area after alcohol application.*



ATR-9N

ATR-9N Wide Cart Trellis Rail System

ATR-9X Additional Wide Cart Hanging Bar

ANS-9 Narrow Shelf for Wide Cart

Dim.: 25.4cm x 87.6cm. Adds up to 53.3cm of height to cart.
Wt. capacity: 23kg



ANS-9



ASH-4

ScopeSafe™

Key-Locking Storage for Delicate Scopes

ASH-4 ScopeSafe™ Locking Scope Holder

- Holds four scopes
- Has a padded interior for cushioned protection
- Contains multiple vents for ample ventilation
- Contains a protective plexiglas door front
- Attaches easily to wall or cart

- Built with a durable steel frame
 - Overall size: 122cm High x 15.2cm Deep x 40.6cm Wide
- Weight: 20 kg.

Takes two Accessory Panels to mount to a Premier Cart. Not compatible with Mini or Narrow Carts.

ASH-2 ScopeSafe™ Locking Scope Holder (not pictured)

- Holds two scopes
 - Same high-quality features as ASH-4
 - Compact: 122cm High x 15.2cm Deep x 22.9cm Wide
- Weight: 14.1 kg.

Takes two Accessory Panels to mount to a Premier Cart. Two ASH-2 ScopeSafes can be mounted to one side of a cart. Not compatible with Mini or Narrow Carts.

Please note: ASH-2 and ASH-4 will not hold colonoscopes.

Isolation Carts



AKL-Y-400IP

AKL-Y-400IP **STANDARD** Steel Key Locking Four-Drawer Isolation Cart Package

Includes AKL-Y-4 Four-Drawer Cart, PSR-1 Tubular Railing, APT-2 Plastic Top, AUH-1 Utility Hooks, and ACR-1 Coat Rack.

For Standard Cart features, see page 4.

For Standard Cart dimensions, see page 6.



PKL-Y-4

PKL-Y-4 Basic **PREMIER™** Aluminum Key Locking Four-Drawer 76cm Isolation Cart

PBL-Y-4 Basic **PREMIER™** Aluminum Breakaway Locking Four-Drawer 76cm Isolation Cart
Drawers: 7.3cm, 23.8cm (3)

For Premier™ Cart features, see page 5.

For Premier™ Cart dimensions, see page 7.



AMC-1Y

AMC-1Y Five-Drawer 61cm Key Locking **STEEL** Mini Cart

Dimensions:

Overall: 88.3cm High x 52.4cm Deep x 54.6cm Wide.

Drawer size: 43.2cm Deep x 38.1cm Wide.

Drawer heights: 7.3cm (3), 15.6cm, 23.8cm.

Weight: 47kg



AMC-2Y

AMC-2Y Six-Drawer 76cm Key Locking **ALUMINUM** Narrow Cart

Dimensions:

Overall: 111.8cm High x 62.9cm Deep x 62.2cm Wide (includes bumper & casters).

Drawer size: 43.2cm Deep x 38.1cm Wide.

Drawer heights: 7.3cm (3), 15.6cm (2), 23.8cm (1).

Weight: 38.5kg



AKL-Y-4

AKL-Y-4 Basic **STANDARD** Steel Key Locking Four-Drawer 76cm Isolation Cart

Drawers: 7.3cm, 23.8cm (3)

For Standard Cart features, see page 4.

For Standard Cart dimensions, see page 6.



ABL-Y-3

ABL-Y-3 Basic **STANDARD** Steel Breakaway Locking Three-Drawer 61cm Isolation Cart

Drawers: 15.6cm, 23.8cm (2)

For Standard Cart features, see page 4.

For Standard Cart dimensions, see page 6.


Cart Accessories



ACG-1 Caster Guard, 4/set
For 12.7cm dia. wheels. Durable with living hinge, locking tab, closure positioning/support pins.



ATC-3 Tracking Swivel Caster



APH-1 STANDARD Push Handle
Does not fit A-SMART Premier Carts.



ATW-3 Soft-grip Twin Push Handles



ACC-24 Navy Cart Cover for 61cm A-SMART® Cart
ACC-27 Navy Cart Cover for 69cm A-SMART Cart
ACC-30 Navy Cart Cover for 76cm A-SMART Cart
Nylon covers with Velcro tabs.



AQC-1 Silent Knight Tablet Crusher
Includes fifty pouches.
AQC-2 Silent Knight Pouch, 1000/box



APS-5 Plastic Seals, 100/pkg.
Consecutively numbered.



ADL-1 Drawer Liner
Black vinyl and latex free. 43.2cm x 55.9cm. Fits standard size drawers.



ALH-1 Label Holders and Cards (blank), 100/pkg. 7.6cm x 2.5cm
ALH-2 Label Holders, 100/pkg. 17.8cm x 5.1cm
ALH-12 Label Holders, 12/pkg. 17.8cm x 5.1cm



APS-1 Pull-tite Seals, 100/pkg.
25.1cm polypropylene seals open without using scissors. Consecutively numbered.



AUH-1P PREMIER Utility Hooks, 2/pkg.
Mount in corner tracks.



AUH-1 STANDARD Utility Hooks, 2/pkg.
Mount to sides or back of cart.

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see page 34.




ACS-2

***ACS-2 Short Catheter Holder**

- Has metal construction with neutral beige finish
- Requires no screws to mount
- 48.3cm High x 20.3cm Wide x 10.2cm Deep
- Divided into two compartments with non-removable divider

***ACS-3 Short Catheter Holder with Clear Front**

- Clear front enables you to see inventory
- Front panel removes for easy cleaning
- Requires no screws to mount
- 48.3cm High x 10.2cm Wide x 10.2cm Deep



ACS-3



AOB-1

***ACL-2 Long Catheter Holder**

- Fits A-SMART® Carts with 76cm drawer space
- Has metal construction with neutral beige finish
- Requires no screws to mount
- 71.1cm High x 20.3cm Wide x 10.2cm Deep
- Divided into two compartments with non-removable divider

***ACL-3 Long Catheter Holder with Clear Front**

- Fits A-SMART Carts with 76cm drawer space
- Clear front enables you to see inventory
- Front panel removes for easy cleaning
- Requires no screws to mount
- 71.1cm High x 10.2cm Wide x 10.2cm Deep



ACL-2



ADO-1

ADO-1 Divided Organizer with Lid

Dimensions: 8.3cm High x 18.1cm Wide x 18.1cm Deep.

***AOB-1 Oxygen Tank Bracket**

Fits "D" and "E" sized tanks.

AFT-1 STANDARD Fluid Tray System for Cart with 61cm Drawer Space, 3 trays/pkg.

AFT-1P PREMIER™ Fluid Tray System (all Cart sizes), 3 trays/pkg. (not pictured)

AFT-2 STANDARD Fluid Tray System for Cart with 76cm Drawer Space, 3 trays/pkg.

Three height-adjustable trays on rail system. Inner tray dimensions: 43.5cm Long x 11.4cm Deep x 10.2cm High.



AFT-1



PSB-1

***PSB-1 PREMIER™ Sharps Container Brackets**

Fits Sharps containers that are 4.7 liters., 7.6 liters, or 11.3 liters

ASB-1N STANDARD Sharps Container Bracket

Fits Sharps containers that are 4.7 liters., 7.6 liters, or 11.3 liters

**May require PHP-1 panels to mount to a Premier™ Cart.*



ASB-1N

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see page 34.


Cart Accessories



ASB-2

ASB-2 Armstrong Security Box

Lockable box for portable use or inside a 7.3cm drawer. Designed to fit beside an AC-670S Drug Tray. The security box measures 27.9cm Wide x 40.6cm Deep x 6.8cm High.



ATD-2

ATD-2 Tape Dispenser for Trellis Rail

ATD-3 Tape Dispenser for Cart Top

5 tape brackets, each bracket width is 5.7cm (not pictured)

AND-1 Armstrong Lockable Drawer Cover

Strong, secure cover can be fitted to any size A-SMART® Cart drawer (7.3cm, 15.6cm, 23.8cm, or 31.8cm drawer).



AND-1

AWB-1 Wire Basket*

Plastic-coated wire basket. 38.1cm Long x 30.5cm Wide. 10.2cm High in back, 7.6cm High in front. Sloping top, flat bottom.

AWB-3 Wire Chart Holder*

Plastic-coated wire chart holder. 40.6cm Long x 10.2cm Deep x 25.4cm High.



AWB-1

AWB-5 Sloping Wire Basket*

Plastic-coated wire basket measures 30.5cm square. 20.3cm High in back, 10.2cm High in front. Sloping top, flat bottom.



AWB-3

AWC-1 Waste Container*

12.9 liters. 28.9cm x 21cm x 30.8cm.

AWC-2 Fire Retardant Waste Container*

13.2 liters. 31.1cm x 21cm x 29.2cm.

AWC-3 Large Waste Container*

26.5 liters. 36.5cm x 26cm x 38.1cm.



AWB-5

ACV-1PB Push-Button Locking ClearView Security Drug Box with Key Override

CVS-1 ClearView Shelf

PHP-1 PREMIER™ Accessory Panels, 2/pkg.

40.6cm x 10.2cm.

PHP-5 PREMIER™ Accessory Panels for Back of Cart, 2/pkg.

58.4cm x 10.2cm.

AAB-1 STANDARD Accessory Bar Bracket

Required for hanging some accessories on this page. Attaches to two locations on either side of cart. Size: 45.7cm x 2.5cm. (For A-SMART Premier Carts, substitute PHP-1 Accessory Panel).

PHP-1S PREMIER™ Panel for Cart-mounted new-DUET (not pictured)

AHP-3 STANDARD Panel for Cart-mounted new-DUET (not pictured)

*Requires AAB-1 or PHP-1 to hang (see above)



ACV-1PB with CVS-1



PHP-1



AWC-1



AAB-1

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see this page.



PSR-1


PSR-1 Tubular Railing

Made of electro-polished stainless steel. Extends 8.3cm above the top of the cart.

APT-2 **STANDARD** Plastic Top

Seamless ABS plastic protects cart metal and provides easy-to-clean work surface. Rounded corners mean infection has nowhere to hide.

Dimensions: 47cm x 68.6cm.



APT-2

APT-5 **PREMIER™** Plastic Top (not pictured)

Dimensions: 47cm x 68.6cm.

PIV-1 I.V. Pole with Brackets

Adjustable: extends 24.1cm to 100.3cm above cart top.



PIV-1

ASN-3 **STANDARD** Side Shelf

Simple One-Hand Slow Shelf Retraction

Adds 1.8 square feet of temporary work surface when extended and folds neatly away when not in use.

Dimensions: 45.1 cm Deep x 37.8 cm Width.

Weight Limit: 11.4 kg.



ASN-3

PSN-3 **PREMIER™** Side Shelf,

Simple One-Hand Slow Shelf Retraction

Same as ASN-3 (see above).



ACB-1

ACB-1 **STANDARD** Plastic Cardiac Board with Brackets

Heavy-duty plastic board with handle that mounts to back of cart for easy access.

Dimensions: 39.4cm x 54.6cm x 1.9cm.



PCB-1

PCB-1 **PREMIER™** Cardiac Board with Brackets

Same as ACB-1 (see above).

ASU-1 Suction Unit Holder

Dimensions: 25.4cm High x 48.3cm Wide x 19.7cm Deep.



PSN-3



ASU-1

AGB-1 Glove Box

Mounts to twenty different positions on a Standard Cart. Fasteners included. Beige. Gloves not included.

Dimensions: 25.4cm Long x 14.6cm Wide x 9.5cm Deep.

APC-1 Armstrong Wire Punch Card Tray, White

Will hold up to thirty punch cards per tray.

Dimensions: 17.8cm High x 35.6cm Deep x 17.8cm

Wide. Great for fluid transfer.



APC-1



AGB-1

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see page 34.

Shelving Units



ACD-1N PCD-1

ACD-1N STANDARD Defib Corner Shelf Unit
Same features as PCD-1 (see below).

PCD-1 PREMIER™ Defib Corner Shelf Unit
Adds 24.1cm to height of cart. Swivel top rotates 360° – locks every 45°. Adjustable width: 33-46cm. Depth: 30.5cm. Strap included. Capacity: 11.3kg



PDS-1

PNT-2 PREMIER™ Two-shelf Unit, Narrow,
(not pictured)

Dim.: 26cm x 70cm. Adds up to 53.3cm of height. Units include two uprights and two shelves.



ABS-1N

PNS-1 PREMIER™ One-shelf Unit, Narrow
Dim.: 26cm x 70cm. Adds up to 53.3cm of height. Capacity: 25kg

PDS-1 PREMIER™ One-shelf Unit, Deep
Dim.: 37.1cm x 70cm. Adds up to 53.3cm of height. Capacity: 25kg



PNS-1

PDT-2 PREMIER™ Two-shelf Unit, Deep
(not pictured)

Dim.: 37.1cm x 70cm. Adds up to 53.3cm of height. Units include two uprights and two shelves. Capacity: 25kg.



ATN-1N

ABS-1N STANDARD One-shelf Unit, Deep
Adds 50.8cm to height of cart. Outer dimensions: 37.1cm Deep x 70cm Wide. Capacity: 34kg

ABS-2N STANDARD One-shelf Unit, Narrow
Dim.: 26cm x 70cm. Adds up to 50.8cm of height. Capacity: 25kg



ABS-2N

ATN-1N STANDARD Two-shelf Unit, Deep
Outer dimensions: 37.1cm Deep x 69.7cm Wide. Adds 58.4cm to height of cart. Capacity: 25kg

Cart Trellis Systems
See pages 26

ADS-1N STANDARD Armstrong Defib One-shelf Unit with Swivel Top
Outer dimensions: 35.8cm x 68.3cm. Adds 30.5cm to height of cart. Swivel top can rotate 360°. Includes strap to hold monitor. Capacity: 22.7kg

PDS-2 PREMIER™ Armstrong Defib One-shelf Unit with Swivel Top
Same as ADS-1N (see above).



ADS-1N PDS-2

All products are compatible with both Standard & Premier Carts unless otherwise stated.

To mount accessories to A-SMART® Premier™ Carts, extra accessory panels may be required. For panels, see page 34.

Tilt Bins Organizers

Designed for modular assembly, Tilt Bins organize supplies to match your needs. Clear Tilt Bins leave stock visible for quick identification, yet are securely arranged.

Each bin has a slot for labeling and swings open for easy access. Bins are removable for cleaning or refilling.

NOTE: A Trellis Rail System or Shelving Unit (except ACD-1N, PCD-1, PLS-1, and ALS-1) is required to mount Tilt Bins to an A-SMART Cart. See page 26 for Trellis Rail Systems and page 36 for Shelving Units.

ATB-6X Six-Bin Unit

Each bin measures: 11.4cm High x 13.3cm Deep

ATB-5X Five-Bin Unit

Each bin measures: 16.5cm High x 13.3cm Deep

ATB-4X Four-Bin Unit

Each bin measures: 20.6cm High x 16.8cm Deep

ATB-3X Three-Bin Unit

Each bin measures: 24.1cm High x 19.7cm Deep

Each Tilt Bin Organizer has a total width of 57.2cm.

Tilt Bin Dividers


Divide each bin into two sections.

ATB-4D Dividers for ATB-4X,
4/pkg.

ATB-3D Dividers for ATB-3X,
3/pkg.



All products are compatible with both Standard and Premier™ Carts unless otherwise stated.



A combination of bin sizes can be used to accommodate multiple-sized supplies or equipment.

**Need More Organization
For All Your Supplies?
See pages 38-39.**



Lockable Tilt Bins Organizers



NOTE: A Trellis Rail System or Shelving Unit (except ACD-1N, PCD-1, PLS-1, and ALS-1) is required to mount Tilt Bins to an A-SMART Cart. See page 26 for Trellis Rail Systems and page 36 for Shelving Units. **Accessories (including Lockable Tilt Bins) can't be mounted directly above Lockable Tilt Bins.**

Tilt Bin dimensions shown above.


Key Lock Tilt Bins

ATB-5XK Lockable Five-Bin Unit
ATB-4XK Lockable Four-Bin Unit
ATB-3XK Lockable Three-Bin Unit

Breakaway Lock Tilt Bins

ATB-5XL Lockable Five-Bin Unit
ATB-4XL Lockable Four-Bin Unit
ATB-3XL Lockable Three-Bin Unit

A Place For Everything...




AMT-1N

AMT-1N Modular Tray

Full drawer tray includes two vertical rails and eight adjustable dividers with card holder. Dimensions: 54.6cm Wide x 38.7cm Deep x 5.9cm High.

AMT-2N Modular Tray

Same as the AMT-1 (see above), except with an ampule organizer, eleven adjustable dividers, and attached label holders. Dimensions: 54.6cm Wide x 38.7cm Deep x 5.9cm High



AMT-6N

AMT-3N Full Drawer Tray with Three Ampule Holders

Thirty-three dividers, three ampule trays, and four securing clips. Dimensions: 54.6cm Wide x 38.7cm Deep x 5.9cm High



AMT-2N

AMT-6N Modular Drawer Tray


Durable, plastic modular tray fits in a 15.6cm drawer. Includes two vertical rails and six dividers. Dimensions: 54.6cm Wide x 38.7cm Deep x 12.7cm High.



AMT-7N

AMT-7N Modular Tray


Full drawer tray, eighteen ampule dividers, and seven long dividers. Dimensions: 54.6cm Wide x 38.7cm Deep x 5.9cm High



AMT-3N

AMT-8N Modular Tray

Full drawer tray, twenty-nine ampule dividers, and four long dividers. Dimensions: 54.6cm Wide x 38.7cm Deep x 5.9cm High




AMT-8N

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

Mini/Narrow Cart Trays

AMT-4N Full Drawer Tray for Mini Cart

Full drawer tray with nine ampule dividers, ampule tray, rail, and four adjustable dividers. Fits one per drawer. Dimensions: 36.1cm Wide x 40.6cm Deep x 6.4cm High



AMT-4N

AMT-5N Full Drawer Tray for Mini Cart

Full drawer tray with two rails and four adjustable dividers. Fits one per drawer. Dimensions: 36.1cm Wide x 40.6cm Deep x 6.4cm High



AMT-5N

...Everything In Its Place.

ADV-1N Modular Drawer Divider Kit
Subdivides a 7.3cm drawer to better organize your supplies. Includes two 55.9cm rails, five 39.4cm rails, ten 16.5cm dividers, and three "L" clips.



IMMEDIATE, accurate access

FASTER emergency response

EASIER, better inventory control



ADT-1 Crash Cart Drawer Tray with Lid

Control Your Drawer Inventory

Clear acrylic lid with thumb hole can be secured with locking seal.

Inside dim.: 53.3cm Wide x 38.1cm

Deep x 6.3cm High.

Outside dim.: 55.9cm Wide x 40.6cm

Deep x 7cm High.



ADV-10N Modular Divider Kit for ADT-1 Drawer Tray

Fits inside ADT-1 Drawer Tray. Includes two 52.5cm Long grooved plastic rails, four 35.6cm Long rails, six 10.1cm Long dividers, and six 14cm Long dividers.



ADV-6N Modular Drawer Divider Kit for 15.6cm, 23.8cm, and 30.5cm Drawers

12.7cm high modular kit fits 15.6cm, 23.8cm, or 30.5cm drawers. Includes two 55.9cm rails, four 39.4cm rails, and six 16.8cm dividers.



APS-1 Pull Tite™ Seal, 100/pkg.

Easy-to-use Drawer Inventory Control

25.1cm polypropylene seals fit the ADT-1 Drawer Tray. Consecutively numbered.



ADB-1 Security Bag, 20/pkg.

67.3cm x 47cm. Fits over 7.3cm trays (ADT-1, AMT-1, AMT-2, and AMT-3 Trays). Perforated center shows you when it's been used.

ADB-2 Security Bag, 20/pkg.

82.5cm x 54.6cm. Fits over 12.7cm tray.

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

Drawer Trays for Stock Exchange



Drug Trays

AC-670S Divisible Drug Tray, 7.3cm Deep

Sand color. Inside dimensions: 37.8cm Long x 23.5cm Wide. Holds up to eleven short dividers and seven long dividers (*not included*). Drawers hold a maximum of two trays, except for the Mini and Narrow Carts.

AC-675S Short Divider for AC-670S
Max. slots: 11

AC-680S Long Divider for AC-670S
Max. slots: 7

AC-672S Divisible Drug Tray, 12.7cm Deep

Same as the AC-670S (*see above*), but made for deeper drawers. Drawers hold a maximum of two trays.

AC-685S Short Divider for AC-672S
Max. slots: 11

AC-690S Long Divider for AC-672S
Max. slots: 7

ACC-655S Clear Acrylic Tray Cover, 2/pkg.

See-through cover fits AC-670S and AC-672S Divisible Drug Trays. Thumb notch ensures easy access.

AC-665 Drug Tray, Small, 6.4cm Deep (*not pictured*)

Blue color. Inside dimensions: 23.5cm Long x 16.8cm Wide. Drawers hold a maximum of four trays.

AC-675 Long Divider for AC-665
Max. slots: 5

AC-677 Short Divider for AC-665
Max. slots: 7



Armstrong Security Box

ASB-2 Armstrong Security Box

Double-Lock Security

This Lockable Security Drug Box can be portable or fastened to the inside of a 7.3cm drawer. The Security Box is

designed to fit two per drawer, or next to our AC-670S 7.3cm Drug Tray.

Dimensions: 27.9cm Wide x 40.6cm Deep x 6.7cm High.

Plastic Divider Sets

AEF-3 Plastic Divider Set for 7.3cm Drawer

Two 55.9cm Long rails and six 41.3cm Long x 6.7cm High dividers.



AEF-3L Long Plastic Divider Set for 7.3cm Drawer

Four 55.9cm Long rails and two 41.3cm Long x 6.7cm High dividers.

AEF-6 Plastic Divider Set for 15.6cm Drawer

Two 55.9cm Long rails and six 41.3cm Long x 14.3cm High dividers.

AEF-6L Long Plastic Divider Set for 15.6cm Drawer

Four 55.9cm Long rails and two 41.3cm Long x 14.3cm High dividers.

AEF-9 Plastic Divider Set for 23.8cm Drawer

Two 55.9cm Long rails and six 41.3cm Long x 22.2cm High dividers.



AEF-9L Long Plastic Divider Set for 23.8cm Drawer

Four 55.9cm Long rails and two 41.3cm Long x 22.2cm High dividers.

All products are compatible with both Standard and Premier™ Carts unless otherwise stated.

Cart Extender

ACE-1 Armstrong Cart Extender*, Front-open

- Adds 1.5 cu. meters to cart
- Key lock to protect contents
- Includes three adjustable shelves
- Each shelf fits AC-670S or AC-672S drawer trays (one per shelf)
- Safety foot adds stability
- Dimensions: 87.3cm High x 45.7cm Deep x 34.3cm Wide

**Cart Extenders only fit A-Smart Carts with 76cm of drawer space. Does not fit on a Narrow or Mini Cart.*



Looking For Extra Wide Drawers?

See page 30, Aluminum Push-Button WIDE Cart

Glove Dispensers



AGB-1

AGB-1 Cart-mount Glove Box

14.6cm Wide x 25.4cm High x 9.8cm Deep.

AE-5359 Glove Dispenser, 6/case

Wall-mountable glove dispenser in a light beige. Adjustable shelf fits different sized boxes of gloves.

36.8cm x 15.2cm x 11.4cm.

AGB-2 Steel, Baked-on Enamel Glove Box

Holds 2 boxes of gloves. 25.4cm Wide x 27.3cm High x 9.8cm Deep.

Mounting hardware not included

AGB-3 Steel, Baked-on Enamel Glove Box

40.6cm Wide x 25.4cm High x 9.8cm Deep.



AE-5359

AGB-4 Steel, Baked-on Enamel Glove Box

(not shown)

Holds 4 boxes of gloves.

57.2cm Wide x 25.4cm High x 9.8cm Deep.

Mounting hardware not included

AGB-2P Plexiglas Glove Box

25.4cm Wide x 27.3cm High x 9.8cm Deep.

Mounting hardware not included

AGB-3P Plexiglas Glove Box *(not shown)*

Holds 3 boxes of gloves.

40.6cm Wide x 25.4cm High x 9.8cm Deep.

Mounting hardware not included



AGB-2



AGB-3



AGB-2P



AGB-3V

AGB-3V Steel, Baked-on Enamel Glove Box

Holds 3 boxes of gloves vertically. 14cm Wide x 76.2cm High x 10.2cm Deep.

Mounting hardware not included