



The Next Generation of Exam Tables just got better

proglide300⁺
power exam table

Premium comfort cushion
OneTouch WheelBase[®]
New compact design
650 lb capacity
ADA height



This is the table that:
Helps protect your staff from injuries
Helps protect your facility from damage
Helps you protect patients from infection

proglide300⁺

power exam table



TABLE SPECIFICATIONS

- Height:** Min. 18.5" (± .5") 46.9cm (±1.27cm) at Seat Max. 37" at seat (93.9cm)
- Width:** 28" (71.12cm)
- Length:** 57" – 76" with foot section extended (144.7cm - 193cm)
- Legrest:** 19" W x 13" D (48.26cm W x 33cm D)
- Front Storage Drawer:** 19" W x 12" D (48.26cm W x 30.48cm D)
- Rear Storage Drawer:** 12.5" W x 24.5" D (31.75cm W x 62.23cm D)
- Stirrups:** Four Lateral Positions, Adjustable Length
- Dual Electro-Mechanical Actuator System** for Efficient Speed
- Electrical:** 120 VAC, 60HZ
- International Model** 4040-650-301: 220/240V, 50/60Hz
4070-650-301: 220/240V, 50/60Hz

ACCESSORIES

- 042 Multi-Positioning Patient Assist Armrests/Safety Rails Pair (Field Installable)
- 251 Bierhoff knee crutches
- 467 Mounting Plate for Welch Allen Light # 48955 - Patient Right
- 468 Mounting Plate for Welch Allen Light # 48955 - Patient Left
- 822 Stainless Steel Drain Pan with Spout (Factory installed only)
- 870 Drawer Warmer 110V (Factory Installed only)
- 872 Drawer Warmer 220V (International Models (Factory Installed only)
- 900 Lifting bracket kit 920 –Easy installation kit
Fire retardant foam (CAL133) available

KEY FEATURES

- Power Hi-Lo + Power Backrest
- Power Hi-Lo with Air Spring Assisted Manual Backrest
- One-function Foot Control – Power Hi-Lo
- Two-function Foot Control – Power Hi-Lo, Power Back
- 650lb (294.84kg) Patient Weight Capacity
- OneTouch WheelBase® Integrated Central-Locking Caster System
- Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
- Hand Control with OneTouch Patient Positioning®
2 Programmable Positions & Return-to-Home
- Removable Seamless Upholstered Cushion
- Side Storage Drawer
- Large Front Storage Drawer with Comfort Leg Pad
- Drawer Warmer
- Sturdy Steel Extension Leg Pullout with Seamless Ultra-Comfort Leg Cushion
- Stirrups
- Pelvic Tilt (7°)
- Stainless Steel Drain Pan
- Hospital Grade Electric Receptacle with Child Safety Cover (Patient Left)
- 3 Year Limited Warranty

	proglide300 ⁺ power exam table 4040-650-300	proglide300 power exam table 4070-650-300
Power Hi-Lo + Power Backrest	✓	–
Power Hi-Lo with Air Spring Assisted Manual Backrest	–	✓
One-function Foot Control – Power Hi-Lo	–	✓
Two-function Foot Control – Power Hi-Lo, Power Back	✓	–
650lb (294.84kg) Patient Weight Capacity	✓	✓
OneTouch WheelBase® Integrated Central-Locking Caster System	✓	✓
Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)	✓	✓
Hand Control with OneTouch Patient Positioning® 2 Programmable Positions & Return-to-Home	✓	✓
Removable Seamless Upholstered Cushion	✓	✓
Side Storage Drawer	✓	✓
Large Front Storage Drawer with Comfort Leg Pad	✓	✓
Drawer Warmer	✓	✓
Sturdy Steel Extension Leg Pullout with Seamless Ultra-Comfort Leg Cushion	✓	✓
Stirrups	✓	✓
Pelvic Tilt (7°)	✓	✓
Stainless Steel Drain Pan	✓	✓
Hospital Grade Electric Receptacle with Child Safety Cover (Patient Left)	✓	✓
3 Year Limited Warranty	✓	✓

COLORS



Why wouldn't you want a table that moves?

Features that matter - protecting your patients, your staff, and your practice.

OneTouch WheelBase™

Protect your patients from infection by quickly and easily cleaning around and under the exam tables.

Protect your workers and prevent injuries by eliminating the need for heavy lifting and pushing.

Prevent damage to your floors and walls when moving exam tables, and prevent dents in your tables when waxing floors.

Flex your space for your rapidly changing health care needs - Adjust the table placement for different caregiver needs, to make room for a stretcher, or to change the room layout.

OneTouch Patient Positioning™

Featuring a Hand control with two memory positions & Return to Home Feature. Hand and foot controls make it easy for you to focus on your patients.


Additional Features

The 650 lb. (294.84 KG) Patient Weight Capacity and 18.5" (49.5 cm) ADA low access height provides fullest accessibility for your patients.

The standard two large drawers, infinitely adjustable stirrups, dual electrical outlet, and a drawer warmer give you the ability to provide optimal patient-centric care.



OneTouch Locks and Unlocks the base.
Glide your exam table around with ease.



OneTouch WheelBase™



Optional Armrests Rotate to Safety Rail and Transfer Assist Positions

Have you ever seen what's lurking under your exam table?

Infections Cost the Health Care Industry \$10 Billion a Year (cdc.gov)

DO YOU EVER
HAVE BODY
FLUIDS SPILL ON
THE FLOOR?




DO YOU EVER
HAVE SOLUTION
SPILLS OR
DRIPS?

"Vulnerable patient populations rely on frequent and intensive use of outpatient care to maintain or improve their health... It is critical that all of this care be provided under conditions that minimize or eliminate risks of healthcare-associated infections (HAI). All healthcare settings, regardless of the level of care provided, must make infection prevention a priority."

Guide to Infection Prevention For Outpatient Settings: Minimum Expectations for Safe Care, CDC, National Center for Emerging and Zoonotic Diseases. Nov. 2015

The OneTouch WheelBase® is a game changer!

Integrated central-locking wheelbase enables staff to easily move the table for cleaning and lock it solidly into place for exams.



Easy moving
means easy
cleaning

The proGlide is the only Exam Table that can be moved & maintain ADA height.



The **proGlide300+ Exam Table** is designed to optimize the patient and caregiver experience. Because your patient population is changing, every detail matters. That's why UMF Medical is setting a **new benchmark for infection control, smart efficiencies, optimal accessibility, and lasting value.** By asking what's going to improve your patients' experience, UMF Medical has arrived at an exam room solution that is the perfect blend of style and function.

THIS IS THE TABLE YOU'VE BEEN WAITING FOR!

UMFmedical

1316 Eisenhower Boulevard
Johnstown, PA 15904-3393
Phone: 814-266-8726
Toll Free: 800-638-5322
Fax: 814-266-1870

www.UMFmedical.com

