



Head Centric Procedure Chair with Wheelbase®, Programmable Hand and Foot Controls: Model 4010-650-300

SPECIFICATIONS

Height: Min. 18.5" (±.5") 46.9cm (±1.27cm) at seat Max. 37" at seat (93.9cm)
Seat Width: 28" (71.12cm)
Back Width: 24.5" (62.23cm)
Length: 70" – 80" with foot section extended (177.8cm – 203.2cm)
Leg Rest: 28" W x 13" D (71.12cm W x 33.02cm D)
Dual Electro-Mechanical Actuator System for Efficient Speed
Electrical: 120 VAC, 60HZ
International Model 4010-650-301: 220/240V, 50/60Hz
Net Weight: 450lb (204kg)
Shipping Weight: 490lb (222kg)
Shipping Dimensions: L64xW30xH41 (163cm x 76cm x 104cm)

ACCESSORIES

047 Adjustable / Self-Leveling Arms (Pair)
 050 Removable U-Shaped Headrest
 051 Removable Rectangular Headrest
 052 Removable Oval Headrest
 467 Welch Allyn Light Bkt For #48955 (Patient Right)
 468 Welch Allyn Light Bkt For #48955 (Patient Left)
 830 Patient Safety Strap (Factory Installed Only)
 873 220v Outlet Package (Factory Installed) (*International*)
 874 220v Power Cord (*International*)
 900 Lifting Bracket Kit
 910 Easy Power Table Installation Kit
 Fire Retardant Foam (Cal133) Available

KEY FEATURES

- Chair Designed for Head and Neck Centric Practices
- OneTouch Wheelbase®
- 650 lb. (294.84 kg) Patient Weight Capacity
- Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
- 4 Power Functions: Hi-Lo, Back, Tilt, Leg Lift & Extension
- Hand Control with OneTouch Patient Positioning®
2 Programmable Positions & Return-to-Home
- Foot Control with OneTouch Patient Positioning®
2 Programmable Positions and Return-to-Home
- Double-articulating Headrest with Height Adjustability and Interchangeable Headrest Options
- Slim Back Cushion for Accessibility
- Removable Seamless Upholstered Cushion
- Trendelenburg Position
- Hospital-Grade Electrical Receptacle with Child Safety Cover (Patient Left)
- Concealed Paper Roll Holder
- 3 Year Limited Warrant

COLORS

