







Head Centric Procedure Chair with Programmable Hand Control: Model 4010-650-100

SPECIFICATIONS

Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at seat Max. 37" at seat 93.9cm

Seat Width: 28" (71.12cm)

Back Width: 24.5" (62.23cm)

Length: 70" – 80" with foot section extended (177.8cm – 203.2cm)

Leg Rest: 28" W x 13" D (71.12cm W x 33.02cm D)

Dual Electro-Mechanical Actuator System for Efficient Speed Electrical: 120 VAC, 60HZ

International Model 4010-650-101: 220/240V, 50/60Hz

Net Weight: 410lb (185kg)

Shipping Weight: 450lb (204kg)

Shipping Dimensions: L64xW30xH41 (163cm x 76cm x 104cm)

KEY FEATURES

- Chair Designed for Head and Neck Centric Practices
- 650 lb. (294.84 kg) Patient Weight Capacity
- Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
- 4 Power Functions: Hi-Lo, Back, Tilt, Leg Lift & Extension
- Hand Control with OneTouch Patient Positioning[®]
 2 Programmable Positions & Return-to-Home
- Double-articulating Headrest with Height Adjustability and Interchangeable Headrest Options
- Slim Back Cushion for Accessibility
- Removable Seamless Upholstered Cushion
- Trendelenburg Position
- Hospital-Grade Electrical Receptacle with Child Safety Cover (Patient Left)
- Concealed Paper Roll Holder
- 3 Year Limited Warranty

ACCESSORIES

- 047 Adjustable / Self-Leveling Arms (Pair)
- 050 Removable U-Shaped Headrest
- 051 Removable Rectangular Headrest
- 052 Removable Oval Headrest
- 467 Welch Allyn Light Bkt For #48955 (Patient Right)
- 468 Welch Allyn Light Bkt For #48955 (Patient Left)
- 830 Patient Safety Strap (Factory Installed Only)
- 873 220v Outlet Package (Factory Installed) (International)
- 874 220v Power Cord (International)
- 900 Lifting Bracket Kit
- 910 Easy Power Table Installation Kit
- Fire Retardant Foam (Cal133) Available

