







Head Centric Procedure Chair with Programmable Hand and Foot Controls: Model 4010-650-200

SPECIFICATIONS

Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at seat Max. 37" at

seat (93.9cm)

Seat Width: 28" (71.12cm) Back Width: 24.5" (62.23cm)

Length: 70" - 80" with foot section extended (177.8cm -

203.2cm)

Leg Rest: 28" W x 13" D (71.12cm W x 33.02cm D)

Dual Electro-Mechanical Actuator System for Efficient Speed

Electrical: 120 VAC. 60HZ

International Model 4010-650-201: 220/240V, 50/60Hz

Net Weight: 420lb (190kg)

ACCESSORIES

047 Adjustable / Self-Leveling Arms (Pair)

050 Removable U-Shaped Headrest

051 Removable Rectangular Headrest

052 Removable Oval Headrest

467 Welch Allyn Light Bkt For #48955 (Patient Right)

468 Welch Allyn Light Bkt For #48955 (Patient Left)

830 Patient Safety Strap (Factory Installed Only)

873 220v Outlet Package (Factory Installed) (International)

874 220v Power Cord (International)

900 Lifting Bracket Kit

910 Easy Power Table Installation Kit

Fire Retardant Foam (Cal133) Available

KEY FEATURES

- Chair Designed for Head and Neck Centric Practices
- 650 lb. (294.84 kg) Patient Weight Capacity
- Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
- 4 Power Functions: Hi-Lo, Back, Tilt, Leg Lift & Extension
- Hand Control with OneTouch Patient Positioning®
- 2 Programmable Positions & Return-to-Home
- Foot Control with OneTouch Patient Positioning®
- 2 Programmable Positions and Return-to-Home
- Double-articulating Headrest with Height Adjustability and Interchangeable Headrest Options
- Slim Back Cushion for Accessibility
- Removable Seamless Upholstered Cushion
- Trendelenburg Position
- Hospital-Grade Electrical Receptacle with Child Safety Cover (Patient Left)
- Concealed Paper Roll Holder
- 3 Year Limited Warranty

COLORS

















